

## MICROFIFCTROOF HMPHFIFRS with RFCHARGFARIF GRASS

P15 Low noise, differential, capacity compensated <u>AC</u> amplifier for recording nervermuscle action potentials - with spike/EEG filter.

P 17 2 single-sided FET inputs. Unity (0.9) gain DC high impedance probe to interface high impedance source with a differential amplifier or chopper amplifier, i.e., Grass Model 7P1 Preamplifier or CRO. Not capacity compensated.

Low noise, differential, capacity compensated AC-DC amplifier. Guard to drive shield of probe for microelectrode recording of cortical, muscle, nerve action or routine resting potentials. Meter to direct readout of electrode resistance. Simultaneous AC and DC outputs plus many other convenient features.

P18

Low noise, low gain, capacity compensated

DC microelectrode differential amplifier
for routine measurements of nerve-muscle
resting potentials and intracellular potentials.

® GRASS SINCE 1935 MEDICAL INSTRUMENTS

© GRASS INSTRUMENT CO. 1968

P16

D102F8

QUINCY, MA., 02169 U.S.A.